

In Silico – Design and computer validation for PCR probes

This document provides a general checklist of steps to develop a PCR probe, but is not intended to be detailed guidance.

- 1) Obtain a [reference mitochondrial sequence](#) for the target species.
 - a) Prefer full mitochondrial genome, but if that is not available pick a reference sequence of a mitochondrial gene.
 - b) Prefer commonly sequenced regions of common genes, ie COI, 12S... This helps ensure that we are selecting unique sequences, not just unique to the database sequences.
 - c) Prefer genes that we have metabarcoding primers for. This will allow us to use metabarcoding as a directly comparative measure.
- 2) Run [NCBI Primer-BLAST](#), with the reference sequence as the PCR Template
 - a) PCR product size should be set to 70 – 140
 - b) The specificity check should be run against the nr database.
 - i) Organism should be left blank
 - c) If there are multiple primer options which identically match the target, prefer the primer that is:
 - i) On COI or 12S, or if not possible on another commonly used gene.
 - ii) Has the fewest potentially unintended targets.
 - iii) Has the fewest unintended targets closely related / geographically located near to the target species.
 - iv) If there are unavoidable potentially unintended targets, record them for probe design.
- 3) Use [IDT PrimerQuest](#), with the same reference sequence and the primers designed by Primer-BLAST.
 - a) Enter sequence and select qPCR 2 Primers + probe
 - b) Show custom design parameters and enter the chosen forward and reverse primers.
 - c) Run it
 - d) If the forward and reverse primers were not specific to just the target species;
 - i) Take the probe and run a [BLASTn](#) against all of the unintended targets returned by Primer-BLAST
 - ii) If the probe has mismatches with all of the unintended targets you are all set.
 - iii) If the probe does not have mismatches with any of them, use [BLASTn](#) to align the target reference with the troublesome unintended target.
 - (1) In that alignment identify the location of the primer suggested by IDT.
 - (2) Identify the smallest shift of the location of the probe that would introduce a mismatch between the two sequences.
 - (3) Return to “ i) “ with the new shifted probe, make sure the GC content didn’t change too much, and that you didn’t add in something like “GGGG”...

In Vitro – Positive control, calibration, and basic negative control.

- 1) Order a synthetic gene fragment of the region containing the amplified sequence for the target species.
 - a) Use the same reference sequence as was used to design the probes.
 - b) [gBlocks Gene Fragments](#), whatever is the cheapest price that will contain at least the target sequence.
 - c) Select a region of the reference sequence that is the longest possible for the price point, currently 500 bp, try and center the amplified region as much as possible.
 - d) If the ordering tool complains about regions on the edge of the sequence that can be easily removed do so.
- 2) dPCR calibration
 - a) Run positive control, negative control
 - i) Positive control, aim to “saturate” the chip with the synthetic sequence.
 - ii) Negative control, water with no synthetic sequence.
 - b) Run calibration chip
 - i) Dilute synthetic reads to ~10,000 copies for the 20,000 well chip. We want between 10% and 90% occupancy. Amount (ng) = length (bp) * 1.1 * 10⁽⁻⁷⁾
- 3) qPCR calibration based on [Forootan et al.](#)
 - a) Run a positive control, negative control
 - i) Positive control, use “a lot” of DNA
 - ii) Negative control, run water with no synthetic sequence.
 - b) 2 fold dilution series from 1 – 64 molecules, with 24 wells each, except for 1 with 48 wells
 - i) Can be scaled, precision depends on the number of replicates and accuracy depends on the spacing of the dilutions. Forootan et al. suggests the most pragmatic approach is to do a course series to estimate the Limit of Detection, then hone in as needed.

In Vivo – Test with real samples

- 1) Begin with sample(s) taken with as many other methods of detection as well.
 - a) Negative controls
 - i) Lab blank
 - ii) Location blank
 - b) Positive control
 - i) Spike a sample with the reference dna (idea from [Perez et al.](#)).
 - (1) Aim for ~10-20x the Limit of Detection, Perez et al. used 100 molecules.
 - ii) Run / compare to a calibrated pure synthetic positive control to see the inhibition factor.
 - c) Suspected positives, from detection with other methods.
 - i) “bucket” samples where we can be almost positive there is dna in the sample.
 - ii) “Normal” samples
 - (1) what volume is needed to get a positive?
 - (2) How does the sensitivity compare to other sampling methods?
 - d) Suspected negatives, from no detection with other methods.
 - i) Where there is almost certainly none of the species, like a sample from Oregon for an east coast only fish...
 - ii) Where we just did not get any hits from other methods.
 - (1) How does the sensitivity compare to other sampling methods?